

Cardiac Arrest Event Data

At Your Fingertips

Avive's Rapid Incident Reporting enables the quick transfer of valuable emergency data from the Avive Connect AED™ to the people who need it, when they need it.

When EMS arrives on scene, first responders can seamlessly access valuable cardiac arrest data at a summary level from the Avive AED™. In addition, the device also displays an Incident ID which can be used by downstream providers, providing access to detailed data prior to EMS arrival at the hospital. Incident Reports which include ECG recordings, shock data, and usage time stamps, are immediately sent to the cloud post event.



Avive makes accessing valuable incident data simple and fast. No more cumbersome downloads, special cables, CD-ROM software, or additional expense required to access lifesaving data.

“I think a major unmet need in the whole industry has been getting data from AEDs to providers. Getting data from the field is so important and often under-appreciated until you have it and you realize just how valuable it is.”

Dr. Benjamin Abella | Director, Center for Resuscitation Science, UPenn

3 Easy Ways to Access Your Avive AED Incident Data

1 On Your Avive Connect AED

After the emergency, access event data at anytime. Navigate to your **Home screen** > **Your Device** > **Use History** > then select the date to view the Incident ID and QR Code for a rescue attempt.




2 On Your REALConnect™ Platform

Any Administrator or Manager assigned to an AED will have access to its incident data on the REALConnect Platform. Simply login to your REALConnect Platform and search for the event data with the Incident ID.

Go to the **“Incidents”** tab in REALConnect and search for your report with the Incident ID.

Once you click on an incident, you can download the full ECG report by selecting **“View Full Incident Report”** in the top right. Your report will automatically download as a PDF document.



3 Online

Healthcare providers can log in to view a patient’s full Incident Report by visiting avive.life/Event-Data and searching with the Incident ID from any emergency. Easily download a PDF of the Incident Report on the REALConnect Platform and quickly share with any provider.

To learn more about how to retrieve data from your device, visit help.avive.life